



## **STEM-CELLBANKER for Use in Mesenchymal Stem Cells (MSCs) Cryopreservation**

Mesenchymal stem cells (MSCs) are multipotent stem cells that is able to differentiate into a variety of cell types. Differentiations of MSCs into osteoblasts, chondrocytes, adipocytes have been shown in *in vitro* or *in vivo*. Recent study shows MSCs a promising approach for regenerative medicine for a wide range of applications such as application to bone and cartilage repair, regeneration of cardiomyocytes or the like.

A proven record of successful cryopreservation for MSCs using STEM-CELLBANKER has been established in our collaborative study with Graduate School of Comprehensive Human Sciences Tsukuba.

(Publication: Under construction)

### **Materials and Method**

Human umbilical cord blood-derived MSCs was cultured and collected by centrifugation. Cell concentration was adjusted to  $2 \times 10^5$  cells/ml with STEM-CELLBANKER solution. The mixture was transfer to a cryotube and frozen directly in  $-80$  deep freezer for 4 days without programmed freezing, followed by cryopreservation in  $-196$  liquid nitrogen tank for 72 days. After thawing the cells in  $37$  waterbath, cell-washing was conducted in a tube with mixture of DMEM (9ml) and FBS (1ml) for 2 times. Cells were centrifuged 1000rpm/s for 5min at  $4$  and followed by cell counting (for cell viability). Cells with concentration of  $4 \times 10^5$  cells/ml were plated in 35mm dish for cell proliferation determination. Number of living cells was counted at 4<sup>th</sup> day and 8<sup>th</sup> day after culture. At the same time, cells with concentration of  $5 \times 10^4$  were cultured in 24-well plate with 2 differentiation media as following for cell differentiation determination.

	Osteoblast Differentiation Medium	Fat cell Differentiation Medium
Composition of Medium	IMDM + 1v/v% FBS 0.1 $\mu$ M dexamethasone 50 $\mu$ g/ml ascorbic acid 10mM $\beta$ -glycerol phosphate 10ng/ml hEGF	IMDM + 1v/v% FBS 0.1 $\mu$ M dexamethasone 0.5mM 3-isobutyl-1- methylxanthine 0.1mM indomethacin

**Result:**

➤ Cell viability

**An average of more than 90% of cell viability after thawing was observed.**

➤ Cell proliferation ability

Figure 1 shows proliferation of MSCs after thawing.

**Thawed MSCs was indicated for possessing cell proliferation ability as normal.**

➤ Cell differentiation ability

Figure 2a. shows Osteoblast differentiation of MSCs after thawing.

Figure 2b. shows Fat cell differentiation of MSCs after thawing.

**Thawed MSCs was indicated for possessing cell differentiation ability as normal.**

**Figure 1**





Fig. 1 shows the proliferation of MSCs after thawing. Increase of cells number was observed on Day 4 and Day 8 indicating that thawed MSCs possesses cell proliferation ability as normal.

**Figure 2a (Staining with Alizalin Red)**



**Figure 2b (Staining with Oil Red O)**



### **Conclusion**

STEM-CELLBANKER showed as high as 90% cell viability while retaining normal cell proliferation and differentiation ability of MSCs. The procedures for cell freezing and cell thawing are very simple. Cells can be frozen directly in -80 degree Celsius for long-term preservation with high cell viability.

Data provided by **Dr. Masumi NAGANO**,  
Graduate School of Comprehensive Human Sciences and Center for Tsukuba, Advanced  
Research Alliance, University of Tsukuba.