Protocol of STEM-CELLBANKER[®]

A case of cryopreservation of human ES and iPS cells

Freezing

- 1. Examine and make sure the cell culture free of contamination, in healthy situation and proper confluency, etc.
- 2. Transfer cells into cryotube.
- 3. Remove culture medium by pipetting without centrifugation (*1).
- 4. Directly add STEM-CELLBANKER[®] (*2) cooled at 4°C to cells in the cryotube at 500 uL.
- 5. Place the cryotube into a container (*2) for programmable freezing cooled at 4°C.
- 6. Place the container (*3) into -80°C freezer.
- 7. Next day, remove the cryotube from the container and transfer it to -196℃ liquid nitrogen storage tank.

Perform 5. and 6. quickly.

*1. The procedure is also available to collect cells into centrifugation tube and remove culture medium by centrifugation at 300rpm (20G).

*2. Prior to the procedure, Rho-associated kinase (ROCK) inhibitor, Y-27632 (Merck) can be added at 10 uM (final conc.) to improve the viability in some cases.

*3. As programmable freezing container, Mr.Frosty (5100-0001) (Nalgene Nunc) and BICELL DEWAR FLASK (Nihon Freezer) can be used.

Thawing

- 1. Remove the frozen cell from storage and quickly thaw in a 37°C shaking water bath.
- Immediately dilute and gently mix cells with 8~10ml of cell wash and recovery solution CELLOTION[™] (*1) or culture medium.
- 3. Centrifuge the cryotube (*2) and remove the supernatant by aspirator.
- 4. Gently suspend the cells with appropriate volume of complete cell culture medium.
- 5. Continue the further culture procedures according to standard protocols.
- *1. CELLOTIONTM helps to increase the yield of cells after washing cells. http://www.zenoag.jp/cellbanker/ja/cellotion.html

nup.//www.zenoaq.jp/cenbanker/ja/cenotion.nu

*2. 5~10 minutes at 300G recommended.

Storage and Stability:

Store at 2 - 8°C. The product will be stable for 3 years after the date of manufacture at this storage temperature. If storing for 3 months or longer after the bottle opened, separate into aliquot and keep at -20°C. To prevent from deterioration of the product, frequent freeze-thaw method (more than 3 times) should be avoided.